



RNF128 Mouse mAb

Catalog No	YP-mAb-18811
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	RNF128 GRAIL
Protein Name	E3 ubiquitin-protein ligase RNF128 (Gene related to anergy in lymphocytes protein) (RING finger protein 128)
Immunogen	Synthesized peptide derived from human RNF128
Specificity	This antibody detects endogenous levels of RNF128 at Human, Mouse
Formulation	
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	47kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

