



# UBE2J2 Rabbit pAb

<b>Catalog No</b>	YP-Ab-18808
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human,Mouse
<b>Applications</b>	WB
<b>Gene Name</b>	UBE2J2 NCUBE2
<b>Protein Name</b>	Ubiquitin-conjugating enzyme E2 J2 (Non-canonical ubiquitin-conjugating enzyme 2) (NCUBE-2)
<b>Immunogen</b>	Synthesized peptide derived from human UBE2J2
<b>Specificity</b>	This antibody detects endogenous levels of UBE2J2 at Human, Mouse
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Polyclonal, Rabbit,IgG
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	28kD
<b>Cell Pathway</b>	Endoplasmic reticulum membrane ; Single-pass type IV membrane protein .
<b>Tissue Specificity</b>	
<b>Function</b>	Catalyzes the covalent attachment of ubiquitin to other proteins. Seems to function in the selective degradation of misfolded membrane proteins from the endoplasmic reticulum (ERAD).
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

