



PHF15 Rabbit mAb

Catalog No	YP-Ab-18843
Isotype	IgG
Reactivity	Human; Mouse;
Applications	WB
Gene Name	PHF15 JADE2 KIAA0239
Protein Name	Protein Jade-2 (PHD finger protein 15)
Immunogen	Synthesized peptide derived from human PHF15
Specificity	This antibody detects endogenous levels of PHF15 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	87kD
Cell Pathway	
Tissue Specificity	
Function	Scaffold subunit of some HBO1 complexes, which have a histone H4 acetyltransferase activity . Acts as a E3 ubiquitin-protein ligase mediating the ubiquitination and subsequent proteasomal degradation of target protein histone demethylase KDM1A . Also acts as a ubiquitin ligase E3 toward itself. Positive regulator of neurogenesis (By similarity).
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images