



ARHGDIB Rabbit pAb

Catalog No	YP-Ab-18856
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	ARHGDIB GDIA2 GDID4 RAP1GN1
Protein Name	Rho GDP-dissociation inhibitor 2 (Rho GDI 2) (Ly-GDI) (Rho-GDI beta)
Immunogen	Synthesized peptide derived from human ARHGDIB
Specificity	This antibody detects endogenous levels of ARHGDIB at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	22kD
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	Detected in bone marrow, thymus and spleen.
Function	Regulates the GDP/GTP exchange reaction of the Rho proteins by inhibiting the dissociation of GDP from them, and the subsequent binding of GTP to them . Regulates reorganization of the actin cytoskeleton mediated by Rho family members .
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images