




CGI-43 Monoclonal Antibody

Catalog No	YP-mAb-01608
Isotype	IgG
Reactivity	Human;Mouse;Monkey
Applications	WB
Gene Name	CCZ1/C7orf28B
Protein Name	Vacuolar fusion protein CCZ1 homolog/Vacuolar fusion protein CCZ1 homolog-like
Immunogen	The antiserum was produced against synthesized peptide derived from human CG028. AA range:411-460
Specificity	CGI-43 Monoclonal Antibody detects endogenous levels of CGI-43 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CCZ1; C7orf28A; CGI-43; Vacuolar fusion protein CCZ1 homolog; CCZ1B; C7orf28B; Vacuolar fusion protein CCZ1 homolog B; Vacuolar fusion protein CCZ1 homolog-like
Observed Band	56kD
Cell Pathway	Lysosome membrane .
Tissue Specificity	Brain,Skin,Testis,Uterus,
Function	similarity:Belongs to the UPF0550 family.,
Background	similarity:Belongs to the UPF0550 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CGI-43 Monoclonal Antibody