




FoxD4 Monoclonal Antibody

Catalog No	YP-mAb-01721
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	FOXD4
Protein Name	Forkhead box protein D4
Immunogen	Synthesized peptide derived from the Internal region of human FoxD4.
Specificity	FoxD4 Monoclonal Antibody detects endogenous levels of FoxD4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FOXD4; FKHL9; FOXD4A; FREAC5; Forkhead box protein D4; Forkhead-related protein FKHL9; Forkhead-related transcription factor 5; FREAC-5; Myeloid factor-alpha
Observed Band	47kD
Cell Pathway	Nucleus.
Tissue Specificity	Chondrosarcoma,Heart,PCR rescued clones,
Function	similarity:Contains 1 fork-head DNA-binding domain.,
Background	This gene encodes a member of the forkhead/winged helix-box (FOX) family of transcription factors. FOX transcription factors play critical roles in the regulation of multiple processes including metabolism, cell proliferation and gene expression during ontogenesis. Mutations in this gene are associated with a complex phenotype consisting of dilated cardiomyopathy, obsessive-compulsive disorders, and suicidality. [provided by RefSeq, Mar 2012],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using FoxD4 Monoclonal Antibody