



PRPF39 Monoclonal Antibody

Catalog No	YP-mAb-01963
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PRPF39
Protein Name	Pre-mRNA-processing factor 39
Immunogen	The antiserum was produced against synthesized peptide derived from human PRPF39. AA range:263-312
Specificity	PRPF39 Monoclonal Antibody detects endogenous levels of PRPF39 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PRPF39; Pre-mRNA-processing factor 39; PRP39 homolog
Observed Band	74kD
Cell Pathway	Nucleus .
Tissue Specificity	Ileal mucosa,Liver,Placenta,Uterus,
Function	function:Involved in pre-mRNA splicing.,similarity:Belongs to the PRP39 family.,similarity:Contains 7 HAT repeats.,
Background	function:Involved in pre-mRNA splicing.,similarity:Belongs to the PRP39 family.,similarity:Contains 7 HAT repeats.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

