




ZBTB45 Monoclonal Antibody

Catalog No	YP-mAb-03499
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZBTB45
Protein Name	Zinc finger and BTB domain-containing protein 45
Immunogen	Recombinant Protein of ZBTB45
Specificity	The antibody detects endogenous ZBTB45 protein.
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Zinc finger and BTB domain-containing protein 45; Zinc finger protein 499
Observed Band	70kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,Cervix,Mammary gland,
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 4 C2H2-type zinc fingers.,
Background	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 4 C2H2-type zinc fingers.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ZBTB45 Monoclonal Antibody