




FAM3D Monoclonal Antibody

Catalog No	YP-mAb-03867
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	FAM3D
Protein Name	Protein FAM3D
Immunogen	Synthesized peptide derived from the C-terminal region of human FAM3D.
Specificity	FAM3D Monoclonal Antibody detects endogenous levels of FAM3D protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FAM3D; Protein FAM3D
Observed Band	25kD
Cell Pathway	Secreted .
Tissue Specificity	Abundantly expressed in placenta and weakly expressed in small intestine.
Function	similarity:Belongs to the FAM3 family.,tissue specificity:Abundantly expressed in placenta and weakly expressed in small intestine.,
Background	similarity:Belongs to the FAM3 family.,tissue specificity:Abundantly expressed in placenta and weakly expressed in small intestine.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using FAM3D Monoclonal Antibody