




MRP-L11 Monoclonal Antibody

Catalog No	YP-mAb-03974
Isotype	IgG
Reactivity	Human;Monkey
Applications	WB
Gene Name	MRPL11
Protein Name	39S ribosomal protein L11 mitochondrial
Immunogen	The antiserum was produced against synthesized peptide derived from human MRPL11. AA range:21-70
Specificity	MRP-L11 Monoclonal Antibody detects endogenous levels of MRP-L11 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MRPL11; CGI-113; 39S ribosomal protein L11; mitochondrial; L11mt; MRP-L11
Observed Band	21kD
Cell Pathway	Mitochondrion .
Tissue Specificity	Duodenum,Lymph,
Function	similarity:Belongs to the ribosomal protein L11P family.,
Background	This nuclear gene encodes a 39S subunit component of the mitochondrial ribosome. Alternative splicing results in multiple transcript variants. Pseudogenes for this gene are found on chromosomes 5 and 12. [provided by RefSeq, May 2014],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using MRP-L11 Monoclonal Antibody