



# MRP-L22 Monoclonal Antibody

<b>Catalog No</b>	YP-mAb-03984
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse;Rat
<b>Applications</b>	WB
<b>Gene Name</b>	MRPL22
<b>Protein Name</b>	39S ribosomal protein L22 mitochondrial
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human MRPL22. AA range:91-140
<b>Specificity</b>	MRP-L22 Monoclonal Antibody detects endogenous levels of MRP-L22 protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Monoclonal, Mouse,IgG
<b>Purification</b>	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-1:2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	MRPL22; MRPL25; RPML25; HSPC158; 39S ribosomal protein L22; mitochondrial; L22mt; MRP-L22; 39S ribosomal protein L25, mitochondrial; L25mt; MRP-L25
<b>Observed Band</b>	24kD
<b>Cell Pathway</b>	Mitochondrion .
<b>Tissue Specificity</b>	Carcinoma,Retina,Skin,Testis,Umbilical cord blood,
<b>Function</b>	similarity:Belongs to the ribosomal protein L22P family.,
<b>Background</b>	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein that belongs to the L22 ribosomal protein family. A pseudogene corresponding to this



gene is found on chromosome 4q. Two transcript variants encoding different isoforms have been found for this gene. [provide

**matters needing attention**

Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

**Products Images**

Western Blot analysis of various cells using MRP-L22 Monoclonal Antibody

