



MRP-L24 Monoclonal Antibody


Catalog No	YP-mAb-03985
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MRPL24
Protein Name	39S ribosomal protein L24 mitochondrial
Immunogen	The antiserum was produced against synthesized peptide derived from human MRPL24. AA range:101-150
Specificity	MRP-L24 Monoclonal Antibody detects endogenous levels of MRP-L24 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MRPL24; 39S ribosomal protein L24; mitochondrial; L24mt; MRP-L24
Observed Band	32kD
Cell Pathway	Mitochondrion .
Tissue Specificity	Adipose tissue,Bone marrow,Ovary,
Function	similarity:Belongs to the ribosomal protein L24P family.,similarity:Contains 1 KOW domain.,
Background	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein which is more than twice the size of its E.coli counterpart (EcoL24). Sequence analysis identified two transcript variants that encode the same protein. [provided by RefSeq, Jul 2008],

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using MRP-L24 Monoclonal Antibody