



NARG1 Monoclonal Antibody

Catalog No	YP-mAb-04029
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	NAA15
Protein Name	N-alpha-acetyltransferase 15 NatA auxiliary subunit
Immunogen	The antiserum was produced against synthesized peptide derived from human NARG1. AA range:221-270
Specificity	NARG1 Monoclonal Antibody detects endogenous levels of NARG1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NAA15; GA19; NARG1; NATH; TBDN100; N-alpha-acetyltransferase 15; NatA auxiliary subunit; Gastric cancer antigen Ga19; N-terminal acetyltransferase; NMDA receptor-regulated protein 1; Protein tubedown-1; Tbdn100
Observed Band	100kD
Cell Pathway	Cytoplasm. Nucleus. Mainly cytoplasmic, nuclear in some cases. Present in the free cytosolic and cytoskeleton-bound polysomes, but not in the membrane-bound polysomes.
Tissue Specificity	Expressed at high levels in testis and in ocular endothelial cells. Also found in brain (corpus callosum), heart, colon, bone marrow and at lower levels in most adult tissues, including thyroid, liver, pancreas, mammary and salivary glands, lung, ovary, urogenital system and upper gastrointestinal tract. Overexpressed in gastric cancer, in papillary thyroid carcinomas and in a Burkitt lymphoma cell line (Daudi). Specifically suppressed in abnormal proliferating blood vessels in eyes of patients with proliferative diabetic retinopathy.
Function	function:The ARD1A-NARG1 complex displays alpha (N-terminal) acetyltransferase activity that may be important for vascular, hematopoietic and neuronal growth and development. Required to control retinal neovascularization in adult ocular endothelial cells. In complex with G22P1 and XRCC5 (Ku80), up-regulates transcription from the osteocalcin promoter.,PTM:Cleaved by caspases during apoptosis, resulting in a stable 35 kDa fragment.,sequence



caution:Contaminating sequence. Potential poly-A sequence.,similarity:Contains 8 TPR repeats.,subcellular location:Mainly cytoplasmic, nuclear in some cases. Present in the free cytosolic and cytoskeleton-bound polysomes, but not in the membrane-bound polysomes.,subunit:Interacts with ARD1A, G22P1, NAT13 and XRCC5.,tissue specificity:Expressed at high levels in testis and in ocular endothelial cells. Also found in brain (corpus callosum), heart, colon

Background

This gene encodes a protein of unknown function. However, similarity to proteins in yeast and other species suggests that this protein may be an N-acetyltransferase. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using NARG1 Monoclonal Antibody