




# RFPL4A Monoclonal Antibody

<b>Catalog No</b>	YP-mAb-04117
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Rat;Mouse;
<b>Applications</b>	WB
<b>Gene Name</b>	RFPL4A
<b>Protein Name</b>	Ret finger protein-like 4A
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human RFPL4A. AA range:155-204
<b>Specificity</b>	RFPL4A Monoclonal Antibody detects endogenous levels of RFPL4A protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Monoclonal, Mouse,IgG
<b>Purification</b>	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-1:2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	RFPL4A; RFPL4; RNF210; Ret finger protein-like 4A; RING finger protein 210
<b>Observed Band</b>	42kD
<b>Cell Pathway</b>	Cytoplasm . Nucleus .
<b>Tissue Specificity</b>	
<b>Function</b>	
<b>Background</b>	The ret finger protein-like 4 gene, Rfpl4, encodes a putative E3 ubiquitin-protein ligase expressed in adult germ cells.
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images



Western Blot analysis of various cells using RFPL4A Monoclonal Antibody