




# TP8L1 Monoclonal Antibody

<b>Catalog No</b>	YP-mAb-05009
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse
<b>Applications</b>	WB
<b>Gene Name</b>	TNFAIP8L1
<b>Protein Name</b>	Tumor necrosis factor alpha-induced protein 8-like protein 1 (TNF alpha-induced protein 8-like protein 1) (TNFAIP8-like protein 1)
<b>Immunogen</b>	Synthesized peptide derived from human protein . at AA range: 40-120
<b>Specificity</b>	TP8L1 Monoclonal Antibody detects endogenous levels of protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Source</b>	Monoclonal, Mouse,IgG
<b>Purification</b>	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-1:2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	20kD
<b>Cell Pathway</b>	Cytoplasm .
<b>Tissue Specificity</b>	High expression detected in most carcinoma cell lines, especially in cells transformed with virus genomes.
<b>Function</b>	similarity:Belongs to the TNFAIP8 family.,
<b>Background</b>	similarity:Belongs to the TNFAIP8 family.,
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images



Western Blot analysis of various cells using TP8L1 Monoclonal Antibody