




# BARH2 Monoclonal Antibody

<b>Catalog No</b>	YP-mAb-05027
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse;Rat
<b>Applications</b>	WB
<b>Gene Name</b>	BARHL2
<b>Protein Name</b>	BarH-like 2 homeobox protein
<b>Immunogen</b>	Synthesized peptide derived from human protein . at AA range: 140-220
<b>Specificity</b>	BARH2 Monoclonal Antibody detects endogenous levels of protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Source</b>	Monoclonal, Mouse,IgG
<b>Purification</b>	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-1:2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	42kD
<b>Cell Pathway</b>	Nucleus .
<b>Tissue Specificity</b>	
<b>Function</b>	function:Potential regulator of neural basic helix-loop-helix genes.,similarity:Belongs to the BAR homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,
<b>Background</b>	function:Potential regulator of neural basic helix-loop-helix genes.,similarity:Belongs to the BAR homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images



Western Blot analysis of various cells using BARH2  
Monoclonal Antibody