



KBRS1 Monoclonal Antibody

| | |
|--------------------|--|
| Catalog No | YP-mAb-05190 |
| Isotype | IgG |
| Reactivity | Human;Mouse |
| Applications | WB |
| Gene Name | NKIRAS1 KBRAS1 |
| Protein Name | NF-kappa-B inhibitor-interacting Ras-like protein 1 (I-kappa-B-interacting Ras-like protein 1) (Kappa B-Ras protein 1) (KappaB-Ras1) |
| Immunogen | Synthesized peptide derived from human protein . at AA range: 50-130 |
| Specificity | KBRS1 Monoclonal Antibody detects endogenous levels of protein. |
| Formulation | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide. |
| Source | Monoclonal, Mouse,IgG |
| Purification | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution | WB 1:500-1:2000 |
| Concentration | 1 mg/ml |
| Purity | ≥90% |
| Storage Stability | -20°C/1 year |
| Synonyms | |
| Observed Band | 21kD |
| Cell Pathway | Cytoplasm . |
| Tissue Specificity | Widely expressed. |
| Function | domain:In contrast to other members of the Ras family, the members of the KapMAB-Ras subfamily do not contain the conserved Gly and Gln residues in positions 13 and 65, which are replaced by Leu residues, and are therefore similar to the constitutively active forms of oncogenic forms of Ras. This suggests that members of this family are clearly different from other small GTPases proteins.,function:Atypical Ras-like protein that acts as a potent regulator of NF-kappa-B activity by preventing the degradation of NF-kappa-B inhibitor beta (NFKBIB) by most signals, explaining why NFKBIB is more resistant to degradation. May act by blocking phosphorylation of NFKBIB and mediating cytoplasmic retention of p65/RELA NF-kappa-B subunit. It is unclear whether it acts as a GTPase. Both GTP- and GDP-bound forms block phosphorylation of NFKBIB.,similarity:Belongs to the small GTPase superfamily. Ras f |
| Background | domain:In contrast to other members of the Ras family, the members of the KapMAB-Ras subfamily do not contain the conserved Gly and Gln residues in positions 13 and 65, which are replaced by Leu residues, and are therefore similar |



to the constitutively active forms of oncogenic forms of Ras. This suggests that members of this family are clearly different from other small GTPases proteins.,function:Atypical Ras-like protein that acts as a potent regulator of NF-kappa-B activity by preventing the degradation of NF-kappa-B inhibitor beta (NFKBIB) by most signals, explaining why NFKBIB is more resistant to degradation. May act by blocking phosphorylation of NFKBIB and mediating cytoplasmic retention of p65/RELA NF-kappa-B subunit. It is unclear whether it acts as a GTPase. Both GTP- and GDP-bound forms block phosphorylation of NFKBIB.,similarity:Belongs to the small GTPase superfamily. Ras family. KapMAB-Ras subfamily.,subunit:Interacts with both NF-kappa-B inhibitor alpha (NFKBIA) and beta (NFKBIB) in vitro. However, it probably only interacts with NFKBIB in vivo. Forms a complex with NFKBIB and NF-kappa-B heterodimer (p50/NFKB1 and p65/RELA). Also interacts with c-Rel (REL),tissue specificity:Widely expressed.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using KBRS1 Monoclonal Antibody