




ABHDD Monoclonal Antibody

Catalog No	YP-mAb-05263
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ABHD13 C13orf6
Protein Name	Abhydrolase domain-containing protein 13 (EC 3.-.-.-)
Immunogen	Synthesized peptide derived from human protein . at AA range: 120-200
Specificity	ABHDD Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	37kD
Cell Pathway	Membrane ; Single-pass type II membrane protein .
Tissue Specificity	Brain,Placenta,
Function	sequence caution:Translated as Leu.,similarity:Belongs to the serine esterase family.,
Background	sequence caution:Translated as Leu.,similarity:Belongs to the serine esterase family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ABHDD Monoclonal Antibody