




MLF2 Monoclonal Antibody

Catalog No	YP-mAb-05766
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MLF2
Protein Name	Myeloid leukemia factor 2 (Myelodysplasia-myeloid leukemia factor 2)
Immunogen	Synthesized peptide derived from human protein . at AA range: 120-200
Specificity	MLF2 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	27kD
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Ubiquitously expressed.
Function	similarity:Belongs to the MLF family.,tissue specificity:Ubiquitously expressed.,
Background	similarity:Belongs to the MLF family.,tissue specificity:Ubiquitously expressed.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using MLF2 Monoclonal Antibody