




RAB6B Monoclonal Antibody

Catalog No	YP-mAb-06057
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RAB6B
Protein Name	Ras-related protein Rab-6B
Immunogen	Synthesized peptide derived from human protein . at AA range: 80-160
Specificity	RAB6B Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	22kD
Cell Pathway	Golgi apparatus membrane ; Lipid-anchor . Endoplasmic reticulum-Golgi intermediate compartment . Cytoplasmic vesicle . Colocalizes with BICD1 at vesicular structures that align along microtubules. .
Tissue Specificity	Predominantly expressed in brain.
Function	function:Seems to have a role in retrograde membrane traffic at the level of the Golgi complex.,similarity:Belongs to the small GTPase superfamily. Rab family.,subunit:Interacts with RAB6KIFL.,tissue specificity:Predominantly expressed in brain.,
Background	function:Seems to have a role in retrograde membrane traffic at the level of the Golgi complex.,similarity:Belongs to the small GTPase superfamily. Rab family.,subunit:Interacts with RAB6KIFL.,tissue specificity:Predominantly expressed in brain.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using RAB6B Monoclonal Antibody