




TIPRL Monoclonal Antibody

Catalog No	YP-mAb-06282
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB
Gene Name	TIPRL
Protein Name	TIP41-like protein (Putative MAPK-activating protein PM10) (Type 2A-interacting protein) (TIP)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	TIPRL Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	29kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Fetal kidney,Lung fibroblast,Muscle,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the TIP41 family.,
Background	TIPRL is an inhibitory regulator of protein phosphatase-2A (PP2A) (see PPP2CA; MIM 176915), PP4 (see PPP4C; MIM 602035), and PP6 (see PPP6C; MIM 612725) (McConnell et al., 2007 [PubMed 17384681]).[supplied by OMIM, Nov 2010],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using TIPRL
Monoclonal Antibody