




TPC2L Monoclonal Antibody

Catalog No	YP-mAb-06286
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB
Gene Name	TRAPPC2L HSPC126
Protein Name	Trafficking protein particle complex subunit 2-like protein
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	TPC2L Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	15kD
Cell Pathway	Cytoplasm, perinuclear region . Endoplasmic reticulum . Golgi apparatus .
Tissue Specificity	Expressed in testis, liver, bladder, lung, spleen and brain, several cell lines and primary chondrocytes cell line.
Function	similarity:Belongs to the TRAPP small subunits family. Sedlin subfamily.,
Background	This gene encodes a protein that interacts with the tethering factor trafficking protein particle (TRAPP complex). TRAPP complexes mediate the contact between vesicles and target membranes, and thus, are involved in vesicle-mediated transport of proteins and lipids. The encoded protein is related to the X-linked trafficking protein particle complex 2. A related pseudogene is located on the X chromosome. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jan 2016],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using TPC2L Monoclonal Antibody