




IGFL4 Monoclonal Antibody

Catalog No	YP-mAb-06631
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	IGFL4
Protein Name	Insulin growth factor-like family member 4
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	IGFL4 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	13kD
Cell Pathway	Secreted .
Tissue Specificity	Detected in the cerebellum.
Function	similarity:Belongs to the IGFL family.,
Background	IGF like family member 4(IGFL4) Homo sapiens IGFL4 belongs to the insulin-like growth factor (IGF; see MIM 147440) family of signaling molecules that play critical roles in cellular energy metabolism and in growth and development, especially prenatal growth (Ertage et al., 2006 [PubMed 16890402]).[supplied by OMIM, Mar 2008],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using IGFL4
Monoclonal Antibody