




KLF14 Monoclonal Antibody

Catalog No	YP-mAb-06686
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	KLF14 BTEB5
Protein Name	Krueppel-like factor 14 (Basic transcription element-binding protein 5) (BTE-binding protein 5) (Transcription factor BTEB5)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	KLF14 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	35kD
Cell Pathway	Nucleus .
Tissue Specificity	Epithelium,
Function	similarity:Belongs to the Sp1 C2H2-type zinc-finger protein family.,similarity:Contains 3 C2H2-type zinc fingers.,
Background	This intronless gene encodes a member of the Kruppel-like family of transcription factors. The encoded protein functions as a transcriptional co-repressor, and is induced by transforming growth factor-beta (TGF-beta) to repress TGF-beta receptor II gene expression. This gene exhibits imprinted expression from the maternal allele in embryonic and extra-embryonic tissues. [provided by RefSeq, Jul 2013],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using KLF14 Monoclonal Antibody