




OLIG3 Monoclonal Antibody

Catalog No	YP-mAb-06861
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	OLIG3 BHLHB7 BHLHE20
Protein Name	Oligodendrocyte transcription factor 3 (Oligo3) (Class B basic helix-loop-helix protein 7) (bHLHb7) (Class E basic helix-loop-helix protein 20) (bHLHe20)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	OLIG3 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	29kD
Cell Pathway	Nucleus .
Tissue Specificity	Ovary,
Function	function:May determine the distinct specification program of class A neurons in the dorsal part of the spinal cord and suppress specification of class B neurons.,similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,
Background	function:May determine the distinct specification program of class A neurons in the dorsal part of the spinal cord and suppress specification of class B neurons.,similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using OLIG3 Monoclonal Antibody