




PYGO2 Monoclonal Antibody

Catalog No	YP-mAb-07039
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	PYGO2 PP7910
Protein Name	Pygopus homolog 2
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	PYGO2 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	44kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,Epithelium,Ovary,Placenta,
Function	function:Involved in signal transduction through the Wnt pathway.,sequence caution:Sequencing errors.,similarity:Contains 1 PHD-type zinc finger.,subunit:Binds to BCL9 via the PHD-type zinc finger motif, and thereby becomes part of the nuclear beta-catenin/TCF complex.,
Background	function:Involved in signal transduction through the Wnt pathway.,sequence caution:Sequencing errors.,similarity:Contains 1 PHD-type zinc finger.,subunit:Binds to BCL9 via the PHD-type zinc finger motif, and thereby becomes part of the nuclear beta-catenin/TCF complex.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using PYGO2 Monoclonal Antibody