




T4S18 Monoclonal Antibody

Catalog No	YP-mAb-07086
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	TM4SF18
Protein Name	Transmembrane 4 L6 family member 18
Immunogen	Synthesized peptide derived from human protein . at AA range: 20-100
Specificity	T4S18 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	22kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Kidney,
Function	similarity:Belongs to the L6 tetraspanin family.,
Background	similarity:Belongs to the L6 tetraspanin family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using T4S18 Monoclonal Antibody