




SLAF8 Monoclonal Antibody

Catalog No	YP-mAb-07210
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SLAMF8 BLAME
Protein Name	SLAM family member 8 (B-lymphocyte activator macrophage expressed) (BCM-like membrane protein) (CD antigen CD353)
Immunogen	Synthesized peptide derived from human protein . at AA range: 110-190
Specificity	SLAF8 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	31kD
Cell Pathway	Membrane; Single-pass type I membrane protein.
Tissue Specificity	Expressed in lymph node, spleen, thymus and bone marrow.
Function	function:May play a role in B-lineage commitment and/or modulation of signaling through the B-cell receptor.,similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain.,tissue specificity:Expressed in lymph node, spleen, thymus and bone marrow.,
Background	This gene encodes a member of the CD2 family of cell surface proteins involved in lymphocyte activation. These proteins are characterized by Ig domains. This protein is expressed in lymphoid tissues, and studies of a similar protein in mouse suggest that it may function during B cell lineage commitment. The gene is found in a region of chromosome 1 containing many CD2 genes. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using SLAF8 Monoclonal Antibody