




F176A Monoclonal Antibody

Catalog No	YP-mAb-07331
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	FAM176A TMEM166 SP24
Protein Name	Protein FAM176A (Transmembrane protein 166)
Immunogen	Synthesized peptide derived from human protein . at AA range: 91-140
Specificity	F176A Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	16kD
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass membrane protein . Lysosome membrane ; Single-pass membrane protein .
Tissue Specificity	Expressed in lung, kidney, liver, pancreas, placenta, but not in heart and skeletal muscle.
Function	function:Acts as a regulator of programmed cell death, mediating both autophagy and apoptosis.,similarity:Belongs to the FAM176 family.,tissue specificity:Expressed in lung, kidney, liver, pancreas, placenta, but not in heart and skeletal muscle.,
Background	function:Acts as a regulator of programmed cell death, mediating both autophagy and apoptosis.,similarity:Belongs to the FAM176 family.,tissue specificity:Expressed in lung, kidney, liver, pancreas, placenta, but not in heart and skeletal muscle.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using F176A Monoclonal Antibody