




RN185 Monoclonal Antibody

Catalog No	YP-mAb-07343
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RNF185
Protein Name	E3 ubiquitin-protein ligase RNF185 (EC 6.3.2.-) (RING finger protein 185)
Immunogen	Synthesized peptide derived from human protein . at AA range: 121-170
Specificity	RN185 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	21kD
Cell Pathway	Mitochondrion outer membrane ; Multi-pass membrane protein . Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Ubiquitously expressed.
Function	similarity:Contains 1 RING-type zinc finger.,
Background	similarity:Contains 1 RING-type zinc finger.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using RN185 Monoclonal Antibody