



SCG1 Monoclonal Antibody

Catalog No	YP-mAb-07705
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CHGB SCG1
Protein Name	Secretogranin-1 (Chromogranin-B) (CgB) (Secretogranin I) (SgI) [Cleaved into: GAWK peptide; CCB peptide]
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SCG1 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/ 1 year
Observed Band	
Calculated Molecular Weight	110kD
Cell Pathway	Secreted. Neuroendocrine and endocrine secretory granules.
Tissue Specificity	Expressed in the adrenal medulla, and in pheochromocytoma. Not expressed in liver.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

