



SIA8B Monoclonal Antibody


Catalog No	YP-mAb-07764
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ST8SIA2 SIAT8B STX
Protein Name	Alpha-2,8-sialyltransferase 8B (EC 2.4.99.-) (Sialyltransferase 8B) (SIAT8-B) (Sialyltransferase X) (STX) (Sialyltransferase St8Sia II) (ST8SiaII)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SIA8B Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	41kD
Cell Pathway	Golgi apparatus membrane ; Single-pass type II membrane protein .
Tissue Specificity	Highly expressed in fetal brain, kidney and heart and to a much lesser extent in adult heart and thymus.
Function	function:May transfer sialic acid through alpha-2,8-linkages to the alpha-2,3-linked and alpha-2,6-linked sialic acid of N-linked oligosaccharides of glycoproteins and may be involved in PSA (polysialic acid) expression.,online information:GlycoGene database,online information:ST8Sia II,pathway:Protein modification; protein glycosylation.,similarity:Belongs to the glycosyltransferase 29 family.,tissue specificity:Highly expressed in fetal brain, kidney and heart and to a much lesser extent in adult heart and thymus.,
Background	The protein encoded by this gene is a type II membrane protein that is thought to catalyze the transfer of sialic acid from CMP-sialic acid to N-linked oligosaccharides and glycoproteins. The encoded protein may be found in the Golgi apparatus and may be involved in the production of polysialic acid, a modulator of the adhesive properties of neural cell adhesion molecule (NCAM1). This protein is a member of glycosyltransferase family 29. [provided by RefSeq, Jul 2008],

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using SIA8B Monoclonal Antibody