




BRX1 mouse mAb

Catalog No	YP-mAb-07914
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	BRX1 BRX BXDC2
Protein Name	BRX1
Immunogen	Synthesized peptide derived from human BRX1 AA range: 15-65
Specificity	This antibody detects endogenous levels of BRX1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.28% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Ribosome biogenesis protein BRX1 homolog (Brix domain-containing protein 2)
Observed Band	38kD
Cell Pathway	Nucleus, nucleolus.
Tissue Specificity	Cervix,Placenta,
Function	function:Required for biogenesis of the 60S ribosomal subunit.,similarity:Belongs to the BXDC2/BRX1 family.,similarity:Contains 1 Brix domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using BRX1 mouse mAb