




C2D4D mouse mAb

Catalog No	YP-mAb-07915
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C2CD4D FAM148D
Protein Name	C2D4D
Immunogen	Synthesized peptide derived from human C2D4D AA range: 254-304
Specificity	This antibody detects endogenous levels of C2D4D at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.29% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C2 calcium-dependent domain-containing protein 4D
Observed Band	38kD
Cell Pathway	plasma membrane,
Tissue Specificity	Brain,
Function	similarity:Contains 1 C2 domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using C2D4D mouse mAb