




F132B mouse mAb

Catalog No	YP-mAb-07937
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM132B
Protein Name	F132B
Immunogen	Synthesized peptide derived from human F132B AA range: 102-152
Specificity	This antibody detects endogenous levels of F132B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.51% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein FAM132B
Observed Band	38kD
Cell Pathway	Secreted .
Tissue Specificity	Cerebellum,Testis,
Function	similarity:Belongs to the FAM132 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using F132B mouse mAb