




ZN619 mouse mAb

Catalog No	YP-mAb-07967
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF619
Protein Name	ZN619
Immunogen	Synthesized peptide derived from human ZN619 AA range: 192-242
Specificity	This antibody detects endogenous levels of ZN619 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.81% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Zinc finger protein 619
Observed Band	60kD
Cell Pathway	Nucleus .
Tissue Specificity	Liver,Thyroid,Trachea,Uterus,
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 10 C2H2-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ZN619 mouse mAb