




MET16 mouse mAb

Catalog No	YP-mAb-07982
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	METTL16 METT10D
Protein Name	MET16
Immunogen	Synthesized peptide derived from human MET16 AA range: 449-499
Specificity	This antibody detects endogenous levels of MET16 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.96% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Methyltransferase-like protein 16 (EC 2.1.1.-) (Methyltransferase 10 domain-containing protein)
Observed Band	60kD
Cell Pathway	Nucleus . Cytoplasm .
Tissue Specificity	Muscle,Spinal cord,Uterus,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence caution:Translated as Lys.,similarity:Belongs to the methyltransferase superfamily. METT10D/rImF family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using MET16 mouse mAb