




ZN211 mouse mAb

Catalog No	YP-mAb-08009
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF211
Protein Name	ZN211
Immunogen	Synthesized peptide derived from human ZN211 AA range: 229-279
Specificity	This antibody detects endogenous levels of ZN211 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.123% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Zinc finger protein 211 (Zinc finger protein C2H2-25)
Observed Band	60kD
Cell Pathway	Nucleus .
Tissue Specificity	Blood,Brain,Chondrosarcoma,Endothelial cell,Hippocampus,Skin,Tongue,Uterus,
Function	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,function:May be involved in transcriptional regulation.,sequence caution:Contaminating sequence. Potential poly-A sequence.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 12 C2H2-type zinc fingers.,
Background	This gene encodes a protein containing a Kruppel-associated box domain and multiple zinc finger domains. This protein may play a role in developmental processes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2014],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ZN211 mouse mAb