




GGT3 mouse mAb

Catalog No	YP-mAb-08040
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	GGT3P GGT3
Protein Name	GGT3
Immunogen	Synthesized peptide derived from human GGT3 AA range: 300-350
Specificity	This antibody detects endogenous levels of GGT3 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.154% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Putative gamma-glutamyltranspeptidase 3 (GGT 3) (EC 2.3.2.2) (Gamma-glutamyltransferase 3) (Glutathione hydrolase 3) (EC 3.4.19.13) [Cleaved into: Putative gamma-glutamyltranspeptidase 3 heavy chain; Putative gamma-glutamyltranspeptidase 3 light chain]
Observed Band	60kD
Cell Pathway	Membrane ; Single-pass type II membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using GGT3 mouse mAb