




DCBD2 mouse mAb

Catalog No	YP-mAb-08244
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	DCBLD2 CLCP1 ESDN
Protein Name	DCBD2
Immunogen	Synthesized peptide derived from human DCBD2 AA range: 142-192
Specificity	This antibody detects endogenous levels of DCBD2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.359% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Discoidin, CUB and LCCL domain-containing protein 2 (CUB, LCCL and coagulation factor V/VIII-homology domains protein 1) (Endothelial and smooth muscle cell-derived neuropilin-like protein)
Observed Band	85kD
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Highly expressed in testis, heart, skeletal muscle and also in cultured vascular smooth muscle cells.
Function	induction:Increased in lung cancers during the process of tumor progression.,similarity:Contains 1 CUB domain.,similarity:Contains 1 F5/8 type C domain.,similarity:Contains 1 LCCL domain.,tissue specificity:Highly expressed in testis, heart, skeletal muscle and also in cultured vascular smooth muscle cells.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using DCBD2 mouse mAb