




ACTL9 mouse mAb

Catalog No	YP-mAb-08258
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	ACTL9 ACTL7C HSD21
Protein Name	ACTL9
Immunogen	Synthesized peptide derived from human ACTL9 AA range: 86-136
Specificity	This antibody detects endogenous levels of ACTL9 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasmic vesicle, secretory vesicle, acrosome . Cytoplasm, cytoskeleton, perinuclear theca . Localizes predominantly in the equatorial segment of the sperm head and neck regions, with some localization in the acrosomal segment of the head. Colocalizes in the acrosomal and equatorial segments of sperm with ACTL7A. Colocalizes with PLCZ1 in the equatorial segment of the head of capacitated sperm. .
Tissue Specificity	Testis-specific.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ACTL9 mouse mAb