




F216A mouse mAb

Catalog No	YP-mAb-08269
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	FAM216A C12orf24
Protein Name	F216A
Immunogen	Synthesized peptide derived from human F216A AA range: 171-221
Specificity	This antibody detects endogenous levels of F216A at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	caution:It is uncertain whether Met-1 or Met-18 is the initiator.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using F216A mouse mAb