




ARMC1 mouse mAb

Catalog No	YP-mAb-08272
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ARMC1 ARCP
Protein Name	ARMC1
Immunogen	Synthesized peptide derived from human ARMC1 AA range: 203-253
Specificity	This antibody detects endogenous levels of ARMC1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm . Mitochondrion . Mitochondrion outer membrane . Associates with the outer mitochondrion membrane, most likely through its C-terminus (PubMed:31644573). Not integrated into the mitochondrial outer membrane (PubMed:31644573). .
Tissue Specificity	
Function	similarity:Contains 1 ARM repeat.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using ARMC1 mouse mAb