




ASPX mouse mAb

Catalog No	YP-mAb-08287
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ACRV1
Protein Name	ASPX
Immunogen	Synthesized peptide derived from human ASPX AA range: 180-230
Specificity	This antibody detects endogenous levels of ASPX at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasmic vesicle, secretory vesicle, acrosome. Nascent acrosomal vesicle of Golgi phase spermatids.
Tissue Specificity	Testis.
Function	alternative products:Additional isoforms seem to exist,subcellular location:Nascent acrosomal vesicle of Golgi phase spermatids.,tissue specificity:Testis.,
Background	This gene encodes a testis-specific, differentiation antigen, acrosomal vesicle protein 1, that arises within the acrosomal vesicle during spermatogenesis, and is associated with the acrosomal membranes and matrix of mature sperm. The acrosomal vesicle protein 1 may be involved in sperm-zona binding or penetration. Alternatively spliced transcript variants have been described. [provided by RefSeq, Nov 2010],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ASPX mouse mAb