



DCA16 mouse mAb

Catalog No	YP-mAb-08289
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	DCAF16 C4orf30
Protein Name	DCA16
Immunogen	Synthesized peptide derived from human DCA16 AA range: 58-108
Specificity	This antibody detects endogenous levels of DCA16 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using DCA16 mouse mAb

