




CJ088 mouse mAb

Catalog No	YP-mAb-08292
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	C10orf88
Protein Name	CJ088
Immunogen	Synthesized peptide derived from human CJ088 AA range: 129-179
Specificity	This antibody detects endogenous levels of CJ088 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm . Mitochondrion .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using CJ088 mouse mAb