




BIRC8 mouse mAb

Catalog No	YP-mAb-08307
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	BIRC8 ILP2
Protein Name	BIRC8
Immunogen	Synthesized peptide derived from human BIRC8 AA range: 122-172
Specificity	This antibody detects endogenous levels of BIRC8 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm .
Tissue Specificity	Testis specific in normal tissues.
Function	function:Protects against apoptosis mediated by BAX.,similarity:Belongs to the IAP family.,similarity:Contains 1 BIR repeat.,similarity:Contains 1 RING-type zinc finger.,subunit:Binds to caspase-9.,tissue specificity:Testis specific in normal tissues.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using BIRC8 mouse mAb