




T132D mouse mAb

Catalog No	YP-mAb-08325
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	TMEM132D HBE120 KIAA1944 MOLT
Protein Name	T132D
Immunogen	Synthesized peptide derived from human T132D AA range: 509-559
Specificity	This antibody detects endogenous levels of T132D at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Expressed in mature oligodendrocytes. Detected in the brain, lung, pancreas and testis.
Function	function:May serve as a cell-surface marker for oligodendrocyte differentiation.,similarity:Belongs to the TMEM132 family.,tissue specificity:Expressed in mature oligodendrocytes. Detected in the brain, lung, pancreas and testis.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using T132D mouse mAb