




TM63B mouse mAb

Catalog No	YP-mAb-08352
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TMEM63B C6orf110
Protein Name	TM63B
Immunogen	Synthesized peptide derived from human TM63B AA range: 627-677
Specificity	This antibody detects endogenous levels of TM63B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	similarity:Belongs to the SPO75/TMEM63 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using TM63B mouse mAb