



FAH2B mouse mAb

Catalog No	YP-mAb-08385
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	FAHD2B
Protein Name	FAH2B
Immunogen	Synthesized peptide derived from human FAH2B AA range: 156-206
Specificity	This antibody detects endogenous levels of FAH2B at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	cofactor:Calcium.,cofactor:Magnesium.,function:May have hydrolase activity.,similarity:Belongs to the FAH family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of various cells using FAH2B mouse mAb

